

## REPORT DOCUMENTATION PAGE

AFRL-SR-AR-TR-03-

0370

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, sending the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing the burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Project, Washington, DC 20503.

Reviewing  
Information

1. AGENCY USE ONLY (Leave blank)

2. REPORT DATE

3. REPORT TYPE AND DATES COVERED

01 May 2002 - 30 Apr 2003 Final Report

4. TITLE AND SUBTITLE

DEEP LEVEL TRANSIENT SPECTROSCOPY STUDY OF DEFECTS AND  
IMPURITIES IN GAN AND RELATED III-N MATERIALS

5. FUNDING NUMBERS

61103D

3484/US

6. AUTHOR(S)

DR LOOK

7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)

WRIGHT STATE UNIVERSITY  
3640 COL GELENN HWY  
DAYTON OH 454358. PERFORMING ORGANIZATION  
REPORT NUMBER

9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)

AFOSR/NE  
4015 WILSON BLVD  
SUITE 713  
ARLINGTON VA 2220310. SPONSORING/MONITORING  
AGENCY REPORT NUMBER

F49620-02-1-0214

11. SUPPLEMENTARY NOTES

12a. DISTRIBUTION AVAILABILITY STATEMENT

APPROVED FOR PUBLIC RELEASE, DISTRIBUTION UNLIMITED

12b. DISTRIBUTION CODE

## 13. ABSTRACT (Maximum 200 words)

The new system, for the first time, uses powerful Microsoft Windows Software for various options of DLTS measurements and data analyses. However, Phys Tech, which is a contractor of Accent and responsible for writing the software, did not provide us with a detailed manual of the DLTS Software until July of this year (we only received a draft of the Software manual in July). Therefore, we still have to learn how to operate the system in detail.

20031006 073

14. SUBJECT TERMS

15. NUMBER OF PAGES

16. PRICE CODE

17. SECURITY CLASSIFICATION  
OF REPORT

UNCLASSIFIED

18. SECURITY CLASSIFICATION  
OF THIS PAGE

UNCLASSIFIED

19. SECURITY CLASSIFICATION  
OF ABSTRACT

UNCLASSIFIED

20. LIMITATION OF ABSTRACT

UL

Standard Form 298 (Rev. 2-89) (EG)  
Prescribed by ANSI Std. Z39.18  
Designed using Perform Pro, WHS/DIOR, Oct 94

**Wright State University:**  
**Semiconductor**  
**Research Center**

134 Oelman Hall  
3640 Col. Glenn Hwy.  
Dayton OH 45435

## **Final DURIP Report**

to

**Air Force Office of Scientific Research**

Attn: Dr. Gerald L. Witt  
Physics & Electronics Directorate, AFOSR/NE  
801 N. Randolph St., Arlington, VA 22203-1977  
703-696-8571, Voice  
gerald.witt@afosr.af.mil

**Title: Deep Level Transient Spectroscopy Study of Defects  
and Impurities in GaN and Related III-N Materials**

**Grant F49620-02-1-0214**  
**Principle Investigator: David C. Look**  
**david.look@wpafb.af.mil**

**Date of report: 29Aug03**  
**Period of report: 1May02 – 30April03**  
**Period of grant: 1May02 – 30April03**

**DISTRIBUTION STATEMENT A**  
Approved for Public Release  
Distribution Unlimited

## Report on a new DLTS system (DL8030)

By Zhaoqiang Fang

Semiconductor Research Center, Wright State University, Dayton OH 45435

A new DLTS system (DL8030), which was ordered from Accent Optical Technologies (UK) Ltd., arrived at Wright State University in early October, 2002, and was installed in the second week of December, 2002. The system has two cryostats, i.e., a LN<sub>2</sub> cryostat with four micromanipulators and a closed cycle He cryostat. The DL8030 system is a PC driven, software based, digital Deep Level Transient Spectrometer, which uses the method of Fourier transformation for the investigation of deep level characteristics. We will use the system to study deep levels in wide bandgap semiconductors, i.e., GaN, SiC and ZnO. During the installation and brief training, an engineer from Phys Tech GmbH in Germany showed us many options, which include Capacitance-DLTS, Conductance-DLTS, Isothermal DLTS, and Constant-Capacitance DLTS. The new system, for the first time, uses powerful Microsoft Windows Software for various options of DLTS measurements and data analyses. However, Phys Tech, which is a contractor of Accent and responsible for writing the software, did not provide us with a detailed manual of the DLTS Software until July of this year (we only received a draft of the Software manual in July). Therefore, we still have to learn how to operate the system in detail.

The two examples presented below demonstrate the capability of measuring deep levels below 80 K, by using the He cryostat. We have known that a DLTS center (trap E) at  $E_C - 0.18$  eV can be induced in MOCVD-GaN by 1-MeV electron-irradiation (EI) and the trap actually consists of at least two components, each of which has a thermal energy near 60 meV, but with different capture cross-sections [APL 76, 2086 (2000)]. To fully understand the characteristics of trap E, we have performed 1-MeV EI on an HVPE-GaN sample and measured DLTS spectra (from 300 K to 40 K) using different filling pulse lengths, as shown in Fig. 1. From the figure, we see that the peak of trap E shifts to lower temperatures as the pulse length increases. As compared to traps B and D, which do not shift, trap E has a smaller apparent capture cross-section.

Porous SiC (PSC) has been an object of interest for reducing defects in epitaxial SiC layers grown on PSC substrates. We had earlier performed DLTS measurements on a SiC-on-PSC

sample using our old Bio-Rad DL4600 system; however, many shallow centers below 80 K could not be revealed (Abstracts, 2003 EMC). By using the DL8030 system, Ti and N-related shallow centers (traps E, F, and G) can be clearly observed, as shown in Fig. 2.




Fig. 1 DLTS spectra with different  $W_f$  for HVPE-GaN, irradiated by 1-MeV electrons




Fig. 2 DLTS spectrum for epitaxial 6H-SiC grown on porous 6H-SiC substrate